

MADE IN THE



Milliken  
BioSmart 6

## REUSABLE FACE MASK



COLOR MAY VARY



### ATTENTION

WASH WITH CHLORINE  
BLEACH BEFORE USE  
*see care instructions on back*

P/N: SSP-LBT.RFM

**GENERAL PUBLIC: CAN BE USED UP TO 12 WEEKS ONCE WASHED**

**MEDICAL STAFF: WASH DAILY**

*Effective for up to 75 washes for both the general public and medical staff.*

#### **FOR WASHING:**

##### **MACHINES WITHOUT BLEACH DISPENSER**

1. Use detergent as usual per load requirement.
2. Add 1/2 cup of Clorox or store brand bleach to a normal sized full laundry load. (Or follow directions on label of your bleach product). **\*\*DO NOT USE COLOR SAFE BLEACH or FABRIC SOFTENER.**
3. Add mask/clothes/other items to laundry load (Recommended: To wash with like colors )after adding detergent and bleach.

##### **HIGH EFFICIENCY MACHINES AND OTHER MACHINES WITH DISPENSERS**

1. Use detergent as usual per load requirement.
2. Add Clorox or store brand bleach to fill line in dispenser based on size of the load. (Or follow directions on label of your bleach product). **\*\*DO NOT USE COLOR SAFE BLEACH or FABRIC SOFTENER.**
3. Add mask/clothes/other items to laundry load (Recommended: To wash with like colors ).

***\*\*Bleach activates the fabric properties to aid in blocking the air born virus and bacteria. This material retains the chlorine from the bleach to kill 99% of bacteria and viruses while you are wearing it.***

#### **FOR DRYING:**

1. Air dry or tumble dry on low.

## WEAR INSTRUCTIONS



1) Locate the formable wire intended to mold to your nose. This indicates the top of the mask. Pull the bottom mask elastic over your head and rest on your neck.



2) Raise the bottom elastic on your neck, it should be positioned just below the ears. Locate the top elastic band.



3) Pull the band over your head so it is seated on the crest of your head. If positioned correctly the band should not come into contact with your ears.



4) Place your fingertips on either side of the nose at the top of the mask. Slide your fingertips down both sides of the formable wire to mold the area to the shape of your nose.